








Product Overview (Applications)

Force plate							
Application – Type	9285BA	9260AA	9286BA...	9281E...	9287C...	9253B...	9290BD
Take-off force in performance diagnostics	○	○	○	+	+	+	++
Mobile force measurement (field trials)	—	++	++	+	○	—	+
Highly dynamic force measurement in competition	—	—	○	++	++	+	—
Force measurement under difficult conditions (high load, point load, moisture, etc.)	—	—	+	+	+	++	—
Force measurement under very difficult conditions (highly dynamic, large/small forces, wide measuring range, etc.)	—	—	+	++	++	++	—
Gait and running analysis (varied and dynamic)	+	++	++	++	+	+	—
Gait and balance analysis with very accurate determination of COP, fixed installation	++	++	++	+	+	+	—
Clinical gait and balance analysis with very accurate determination of COP, portable, low installation costs	—	++	++	○	—	—	—